



Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue SARS Antigen Test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current Package Insert instructions and/or Procedure Card included in the kit.

---

**1. How often does the manufacturer recommend that external controls be run on the kit?**

- a) Once every 20 tests
- b) Once per each new shipment of kits
- c) Every first sample tested each day
- d) Once per each untrained operator, once for each new shipment of kits – provided that each different lot in a shipment is tested

---

**2. At what temperature do the kits need to be stored?**

- a) 2°C to 8°C (refrigerated)
- b) 15°C to 30°C (room temperature)
- c) 2°C to 30°C (refrigerated or room temperature)
- d) -20°C (frozen)

---

**3. For proper test performance, which swab should be used when collecting a specimen?**

- a) The swabs provided in the kit.
- b) Any type of swabs
- c) Any Dacron® swabs
- d) Cotton-tipped swabs

---

**4. What different sample types are appropriate for testing with this kit?**

- a) Anterior nares swab only
- b) Anterior nares swab, nasopharyngeal swab, nasal aspirate and nasal wash
- c) Anterior nares swab and nasopharyngeal swab
- d) Anterior nares swab, nasal aspirate and nasal wash

---

**5. How long are samples stable after collection?**

- a) Samples should be tested as soon as possible after collection
- b) Samples can be held in a dry transport tube for up to 5 days at room temperature or 2°C to 8°C.
- b) Both A and B
- b) None of the above

---

**6. How long after adding Reagent Solution into the Reagent Tube is the mixture still good?**

- a) 6 hours
- b) 1 minute
- c) The mixture should be used immediately
- d) 24 hours

---

**7. When should the foil pouch containing the Test Strip be opened?**

- a) At any time
- b) Immediately prior to use
- c) Up to one hour before use
- d) Up to one day before use

---

**8. When should the results on the Test Strip be read?**

- a) 10 minutes
- b) 15 minutes
- c) 5 minutes
- d) 20 minutes

---

**9. Reading the Test Strip before or after the specified read time can produce an inaccurate result.**

- a) True
- b) False

---

**10. If the test result is positive for SARS Antigen, which of the following lines are visible on the Test Strip?**

- a) Only a blue control line
- b) A pink test line and blue control line
- c) Only a pink test line
- d) No lines are visible

Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue SARS Antigen Test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current Package Insert instructions and/or Procedure Card included in the kit.

**1. How often does the manufacturer recommend that external controls be run on the kit?**

- a) Once every 20 tests  
b) Once per each new shipment of kits  
c) Every first sample tested each day  
d) Once per each untrained operator, once for each new shipment of kits – provided that each different lot in a shipment is tested

**2. At what temperature do the kits need to be stored?**

- a) 2°C to 8°C (refrigerated)  
b) 15°C to 30°C (room temperature)  
c) 2°C to 30°C (refrigerated or room temperature)  
d) -20°C (frozen)

**3. For proper test performance, which swab should be used when collecting a specimen?**

- a) The swabs provided in the kit.  
b) Any type of swabs  
c) Any Dacron® swabs  
d) Cotton-tipped swabs

**4. What different sample types are appropriate for testing with this kit?**

- a) Anterior nares swab only  
b) Anterior nares I swab, nasopharyngeal swab, nasal aspirate and nasal wash  
c) Anterior nares swab and nasopharyngeal swab  
d) Anterior nares swab, nasal aspirate and nasal wash

**5. How long are samples stable after collection?**

- a) Samples should be tested as soon as possible after collection  
b) Both A and B  
c) Samples can be held in a dry transport tube for up to 5 days at room temperature or 2°C to 8°C.  
d) None of the above

**6. How long after adding Reagent Solution into the Reagent Tube is the mixture still good?**

- a) 6 hours  
b) 1 minute  
c) The mixture should be used immediately  
d) 24 hours

**7. When should the foil pouch containing the Test Strip be opened?**

- a) At any time  
b) Immediately prior to use  
c) Up to one hour before use  
d) Up to one day before use

**8. When should the results on the Test Strip be read?**

- a) 10 minutes  
b) 15 minutes  
c) 5 minutes  
d) 20 minutes

**9. Reading the Test Strip before or after the specified read time can produce an inaccurate result.**

- a) True  
b) False

**10. If the test result is positive for SARS Antigen, which of the following lines are visible on the Test Strip?**

- a) Only a blue control line  
b) A pink test line and blue control line  
c) Only a pink test line  
d) No lines are visible