Introduction

As part of our ongoing effort to enhance our products, Siemens Healthcare Diagnostics is releasing new software for the DCA Vantage® Analyzer. This customer bulletin provides the upgrade instructions for the new software release.

Once you upgrade to 4.7 you will not be able to downgrade the software. You must have version 4.6 installed before upgrading to 4.7. **If your system is not at version 4.6, please contact your customer support or distributor for more information.**

Upgrading the DCA Vantage Software

Users can download the system software from the Siemens Healthineers Document Library onto a blank USB stick. You then use the USB stick to install software upgrades onto your system.

Upgrading the software takes approximately 10-15 minutes.

The version 4.7 upgrade is only for installation on devices with serial numbers S040000 and higher. Please contact your customer support or distributor if you need help determining your instruments serial number. Users can download system software from the Siemens Healthineers Document Library onto a blank USB stick.

Preparing to Install the Software

Before you can install the software, you need to format a USB stick and then download the software to the USB stick.

You'll need a blank 2, 4, 8, or 16 GB stick formatted to FAT32. Performance testing of this process was performed on SanDisk branded USBs.

You must be a registered user to log in and access the Document Library. If you don't have a login and password, use the link below and follow the instructions before downloading system software:

<https://www.siemens-healthineers.com/services/laboratory-diagnostics/service-and-support/technical-documentation>

Formatting the USB Stick

To format your USB stick, perform the following steps:

1. Insert the USB stick into your computer.
2. Open File Explorer.
3. Right-click on your USB stick.
4. Select **Format**.
5. Open the **File System** dropdown menu and select **FAT32**.

NOTE: Make sure Allocation unit size is set to default.

6. Click **Start**.
7. Once done, click close to go back to your File Explorer folder.

Downloading the Software to the USB Stick

To download the software to the USB stick, perform these steps:

1. Open a web browser on your computer.
2. In the browser go to: <https://doclib.siemens-healthineers.com/home>
3. Log on to the Document Library using your email address and password.
4. In the Document Library Search bar, do one of the following:
 - To search for a specific system, type the name of the system (for example, DCA Vantage 4.7) and click **Enter**.
 - To search for all system updates, type **Software Updates** and click **Enter**.

A new page opens, displaying a list of documentation and software upgrade files matching your search entry.

Note: Search results display in English or your profile default language setting.

5. In the search results list, click the checkbox for the software update you want to download.
6. Insert the blank formatted USB stick into a USB slot on your computer.
7. Click the **Download** button to begin downloading the software upgrade to the download folder of your system. Be sure to record the location of the ZIP file.
8. Unzip the download and copy the files directly onto the formatted USB stick.
9. Continue to the next section to upgrade the software.

Note: Contact your Siemens Healthineers support representative if you have any issues with the download or upgrade installation.

Upgrading the Software

NOTE: You will be prompted for an upgrade code. Before you begin this procedure, contact your technical support provider or distributor.

To start the software upgrade, perform these steps:

1. On the **Home** screen, select **Menu**.
2. Select **System Settings > Additional Settings**.

NOTE: If security is set on the analyzer, to display the **Additional Settings** screen, enter an authorized Supervisor ID or access code.

3. Print the Full System Settings Report. Be sure the printer is enabled and select **Print**.
4. Select **General**.
5. Select the down arrow three times.
6. Select **Install Software**.
7. Insert the USB storage device into the USB port.

CAUTION

Do not remove the USB device from the USB port or turn the analyzer off while the software is installing. These actions can corrupt the data on the analyzer.

8. Select **Next**.

On the **Install Software** screen, a message indicates the system is installing the software.

When the installation completes, a message instructs you to leave the USB device in the USB port and restart the system.

9. To restart the system, select **Restart**.

Finalizing the Software Upgrade

After the DCA Vantage analyzer restarts, if prompted to calibrate the touch screen, follow the onscreen instructions.

1. On the **Install Complete** screen, to display the **Home** screen, select **Back**.
2. Remove the USB device.
3. Print and compare the system settings:
 - a. To print the Full System Settings report, on the **Home** screen, select **Print**.
 - b. Compare the report to the report you printed before the upgrade.
 - c. Change any system settings, if needed.

Software Upgrade Instructions

For more information on the software installation procedures, see the *DCA Vantage Operator's Guide*, Chapter 2 Installation.

If you have questions, contact your technical support provider or distributor.

DCA Vantage is a trademark of Siemens Healthcare Diagnostics.

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