

Durable³, Dependable and Accurate^{1,4} Thermometry



Table of Contents

Thermometry Products

Tympanic Thermometers	3
Genius™ 3 Infrared Tympanic Thermometer	3
Genius [™] Probe Covers and Accessories	4
Electronic Thermometers	6
Filac™ 3000 AD Electronic Thermometer	6
Filac™ 3000 EZ Electronic Thermometer	7
Electronic Thermometer Probe Covers and Accessories	8
Other Electronic Thermometers	10
Other Electronic Thermometers and Probe Covers and Accessories	10
Ultra Forehead Temperature Indicator	11
Miscellaneous Accessories	11

Genius[™] 3 Thermometer: Robust, dependable and accurate¹

- Peak Select technology takes multiple readings and selects the highest temperature for display in 1 to 2 seconds.
- Disposable single-use probe covers help prevent cross-contamination.
- Programmable display shows oral, rectal or an absolute temperature reading in ear mode.
- Redesigned plastic housing allows for an IP22 rating for protection against dust and liquid ingress.⁶

Order Information

Item code	Description	Ship case
303013	Genius™ 3 Thermometer with Base	1

PRODUCT SPECIFICATIONS

Accuracy

• ±0.3°C⁴

Measurement Response Time

• 1–2 seconds

Ambient Operating Range

16°C–33°C (60.8°F–91.4°F)

Temperature Measurement Range

• 33°C-42°C (91.4°F-107.6°F)

Storage Temperature Range

• -25°C−55°C

Battery Life

• 15,000 temperatures (3 AAA batteries)

Anti-Theft System

• Wall mount bracket/Genius[™] mobile cart

Compliant with Standards

- ISO 80601-2-56:2017⁴
- IEC 60601-1:2005/AMD1:20123
- IEC 60601-1-2:2014⁵
- IEC 60601-1-11:2015 Ed.2³

Warranty

• 3 year full manufacturer warranty



Genius™ Tympanic Thermometer Probe Covers

Item code

Description

Disposable, single-use
No-touch disposable probe cover application and ejection
Easy load probe cover
Smooth, rigid plastic with clear membrane
Two probe cover cassettes that fit into base of thermometers
For use with Genius[™] 2 tympanic thermometer

Genius [™] Probe Covers 303030- Cassette of 16 Probe Covers; 6 cassettes/package; 2112 11 packages/box; 2 boxes/case



Ship case

303030



303060

Genius[™] Cart Conversion Kit

• Bracket to attach our Genius[™] 3 thermometer to rolling stand (item #303059)

ltem code	Description	Ship case
303060	Genius™ Cart Conversion Kit	1

Genius[™] Locking Wall Mount

- Mobile thermometer stand with side utility basket or use with Genius[™] 3 tympanic thermometer
- Anti-theft device that helps secure Genius[™] 3 thermometers

Item code	Description	Ship case
303058	Genius™ Locking Wall Mount	1





303097

Genius™ Checker/Calibrator

Calibrates the Genius[™] 2 and Genius[™] 3 tympanic thermometers

Item code	Description	Ship case
303097	Genius™ Checker/Calibrator	1

Filac[™] 3000 AD Electronic Thermometer

- Cold mode for hypothermic patients
- Icon-based user interface or backlit LCD display
- Push button alternates Fahrenheit/Celsius
- Interchangable, color-coded isolation chambers with matching probes
- "Auto on", "Auto off" feature
- Last temperature recall
- 15, 30, 45 and 60 second pulse timer
- Utilizes 202020, 502000 or 500500 disposable probe covers

Order Information

Item code	Description	Ship case
505000	Filac™ 3000 AD Oral/Axillary Complete	1
505005	Filac™ 3000 AD Rectal Complete	1
505027	Filac™ 3000 AD Oral/Axillary Complete – 9 ft	1
505037	Filac™ 3000 AD Rectal Complete – 9 ft	1
500026	Filac™ 3000 Oral Probe – 4 ft	1
500027	Filac™ 3000 Oral Probe – 9 ft	1
500028	Filac™ 3000 Oral Isolation Chamber	1
500036	Filac™ 3000 Rectal Probe – 4 ft	1
500037	Filac™ 3000 Rectal Probe – 9 ft	1
500038	Filac™ 3000 Rectal Isolation Chamber	1
500099	Filac™ 3000 Calibration Plug	1

Filac[™] 3000 AD hardware is warranted for three years from shipment for manufacturer's defect only.





505000

PRODUCT SPECIFICATIONS

Accuracy

- Oral Quick Mode: ±0.3°C
- Oral Standard Mode: ±0.1°C
- Axillary: ±0.1°C
- Rectal: ±0.1°C

Measurement Response Time

- Oral Quick Mode
- 3-5 seconds for non-fever temperatures
 - 8-10 seconds for fever temperatures
- Oral Standard Mode: 6–10 seconds
- Axillary: 8-12 seconds
- Rectal: 10–14 seconds
- Direct Mode: 60-120 seconds

Temperature Measurement Range

• 30°C-43°C (86°F to 109°F)

Battery Life

• 6,000 temperatures (4 AA batteries)

Anti-Theft System

- Electronic anti-theft system
- The Filac[™] 3000 AD Electronic Thermometer comes with the mountable base with key lock (item #500046)

505005

500046

Filac[™] 3000 EZ Electronic Thermometer

- Icon based user interface on backlit LCD display
- Push button alternates Fahrenheit/Celsius
- Interchangeable, color-coded isolation chambers with matching probes
- Auto on", "Auto off" feature
- Last temperature recall
- Utilizes item 202020, 502000, or 500500 disposable probe cover

Order Information

Item code	Description	Ship case
504000	Filac™ 3000 EZ Oral/Axillary Complete	1
504005	Filac™ 3000 EZ Rectal Complete	1
504027	Filac™ 3000 EZ Oral/Axillary Complete – 9ft	1
504037	Filac™ 3000 Rectal Complete – 9 ft	1
500026	Filac™ 3000 Oral Probe – 4 ft	1
500027	Filac™ 3000 Oral Probe – 9ft	1
500028	Filac™ 3000 Oral Isolation Chamber	1
500036	Filac™ 3000 Rectal Probe – 4ft	1
500037	Filac™ 3000 Rectal Probe – 9 ft	1
500038	Filac [™] 3000 Rectal Isolation Chamber	1
500099	Filac [™] 3000 Calibration Plug	1

Filac 3000 EZ hardware is warranted for three years from shipment for manufacturer's defect only.





504005

PRODUCT SPECIFICATIONS

Accuracy

- Oral Quick Mode: ±0.3°C
- Oral Standard Mode: ± 0.1 °C
- Axillary: ±0.1°C
- Rectal: ±0.1°C

Measurement response time

- Oral Quick Mode
- 3-5 seconds for non-fever temperatures
 8-10 seconds for fever temperatures
- 8-10 seconds for lever temperatur
- Axillary: 8–12 seconds
- Rectal: 10–14 seconds
- Direct Mode: 60–120 seconds

Temperature measurement range

• 30-43°C (86°F-109°F)

Battery life

• 6000 temperatures (4 AA batteries)



Filac™ Thermometer Probe Covers

- Disposable, single-use
- Smooth, rigid plastic
- Probe cover package that fits into isolation chamber
- Ejected into waste receptacle by depressing release button
- For use with the Filac[™] 3000 thermometers

ltem code	Description	Ship case
500500	Filac™ Probe Covers 20 probe covers; 25 packages/box; 1 box/case	500
502000	Filac™ Probe Covers 20 probe covers; 25 packages/box; 4 boxes/case	2000
202020	Filac™ Probe Covers 20 probe covers; 25 packages/box; 10 boxes/case	5000



Filac[™] 3000 Thermometer Probes

- Replacement probes for Filac[™] 3000 thermometer
- · Available in oral/axillary or rectal versions
- Available in 4 ft or 9 ft cord lengths

ltem code	Description	Ship case
500026	Filac™ 3000 Oral Probe – 4ft	1
500027	Filac™ 3000 Oral Probe – 9 ft	1
500036	Filac™ 3000 Rectal Probe – 4ft	1
500037	Filac™ 3000 Rectal Probe – 9 ft	1



Filac[™] 3000 Isolation Chamber

- · Available in oral/axillary or rectal versions
- Isolation chamber in house thermometer probes
- Does not include Filac[™] 3000 thermometer probe

ltem code	Description	Ship case
500028	Filac™ 3000 Oral Isolation Chamber	1
500038	Filac™ 3000 Rectal Isolation Chamber	1





500046

Filac[™] 3000 Mounting Base with Pole Clamp or Lock

- Locking base with pole clamp for attachment to mobile stands
- Locking base to secure Filac[™] 3000 thermometer to the wall

ltem code	Description	Ship case
500048	Filac [™] 3000 Mounting Base with Pole Clamp	1
500046	Filac™ 3000 Wall Mounting Base with Lock	1

Filac[™] 3000 Mounting Base for Mobile Stand

- Mounting base to secure Filac[™] 3000 thermometer to mobile stand
- Filac[™] 3000 mounting base with lock (item #500046) not included

ltem code	Description	Ship case
500047	Filac™ 3000 Wall Mounting Base for Mobile Stand	1



Filac[™] 3000 Calibration Plug

• Calibration plug to check calibration of the Filac[™] 3000 thermometer

ltem code	Description	Ship case
500099	Filac [™] 3000 Calibration Plug	1

500099

500047

Other Electronic Thermometers

- Easy-to-read LCD display
- Last temperature taken memory capability
- Audible alarm signals, one beep at startup, 10 beeps max at peak temperature
- Auto turn-off after 5 minutes of inactivity
- Available in Fahrenheit, Celsius or switchable dual scale
- Standard 60 second read or quick read option for faster measurement
- Oral or rectal* thermometers
- Individually boxed or economical polybags available
- Utilizes 16811-PRBO and 16811-PRBR probe covers

*Rectal thermometers come with a red cap for easy identification. Probe covers for rectal thermometers are prelubricated.

Standard Thermometers (~60 second oral read)

Oral Thermometers

ltem code	Description	Ship case
16811-100	Oral, Fahrenheit, Box	144
16811-500	Oral, Fahrenheit, Polybag	500
16811-C	Oral, Celsius, Polybag	144
16811-DS	Oral, Dual Scale, Box	144
16811-EDS	Oral, Dual Scale, Polybag	500
16811-ECO	Oral, Fahrenheit, Polybag	500
Rectal Thermometers		
16811-100R	Rectal, Fahrenheit, Box	72
16811-ECOR	Rectal, Fahrenheit, Polybag	500
16811-500R	Rectal, Fahrenheit, Polybag	500
16811-DSR	Rectal, Dual Scale, Box	72

Quick Read (~10 second oral read)

ltem code	Description	Ship case
16811-EQRO	Quick Read, Oral, Dual Scale, Polybag	500
16811-EQRR	Quick Read, Rectal, Dual Scale, Polybag	500

Probe Covers

Item code	Description	Ship case
16811-PRBO	Oral Probe Covers, Box	2000
16811-PRBR	Rectal Probe Covers, Box	2000



16811-100



Ultra Forehead Temperature Indicator

- · Encapsulated liquid crystal technology, displays accurate temperature changes in seconds
- · Noninvasive and inexpensive alternative to electronic probes
- Indicates temperatures in Celsius and Fahrenheit scales with accuracy to ±1° Fahrenheit
- Ideal for use with laryngeal masks and mac or regional anesthesia when intubation is not used
- · Discrete temperature events: a green color indicates the temperature attained
- Not made with natural rubber latex

ltem code	Description	Ship case
14677-493	Ultra Forehead Temperature Indicator	200

Miscellaneous Accessories



303059

Mobile Stand

- · Portable thermometer stand with probe cover holder
- Includes Genius[™] tympanic thermometer locking mount
- Compatible with the Filac[™] 3000 when also purchased with mounting base for mobile stand (item #500047) and wall mounting base with lock (item #500046)

Item code	Description	Ship case
303059	Mobile Stand with Locking Mount	1

100

102

104

14677-493

- 1. Based on internal report #P-4359-049, revision 3, Ear Thermometer Simulated Use Testing for Accuracy. March 2007.
- 2. Based on internal report #P-4359-049, revision 2, Ear Thermometer Simulated One Year Use Testing for Accuracy. June 2006.
- 3. Based on internal test report #IEC 60601-1, *Medical electrical equipment* Part 1: General requirements for basic safety and essential performance. February 2018.
- 4. Based on internal test report #ISO 80601-2-56, *Medical electrical equipment* Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement. February 2018.
- 5. Based on internal test report #COVI0147, revision 1, Medtronic Genius 3. July 2017.
- 6. Based on internal test report #IEC 60601-1-11, *Medical electrical equipment* Part 1-11: General requirements for basic safety and essential performance. Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment. February 2018.

