



OSOM® Ultra Plus Flu A&B Test

STRONGER CLINICAL PERFORMANCE TAKES LATERAL FLOW TESTING TO THE NEXT LEVEL

Features and Benefits

- Accurate High Sensitivity / Specificity
- CLIA-waived No need for instrument
- Fast and simple workflow Pre-filled Vials, Minimal hands on time, Results in 10 minutes
- Ready-to-use Collection Swabs and Controls Included
- Exceptional Customer Service Online training modules and technical support staffed by certified medical technologists/professionals

Performance Characteristics							
	Sensitivity*	Specificity*	Sample Type	Limit of Detection**	Storage Temperature	Time to Result	
FLU A	90.3%	96.7%	Nasal/NP Swab	7.1x10 ¹ TCID ₅₀ /mL	15-30°C/59-86°F	10 minutes	
FLU B	88.0%	99.2%					








SPECIFICITY





OSOM® Ultra Plus Flu A&B Test

SIMPLE TO USE / RESULTS IN MINUTES



Kit includes one Influenza A+ Control Swab and one Influenza B+ Control Swab

INTERPRETATION OF RESULTS



Ordering Information							
	Configuration	Catalog Number	CPT Codes				
OSOM ULTRA PLUS FLU A&B TEST KIT	25 tests per kit	1032	Flu A: 87804QW				
OSOM ULTRA PLUS FLU A&B CONTROL KIT	5 Flu A Swabs 5 Flu B Swabs	1034	Flu B: 87804QW-59				

To learn more about all of our OSOM Tests scan the code or go to osomtests.com

THE AMERICAS

SEKISUI Diagnostics, LLC One Wall Street Burlington, MA 01803 Phone: 800 332 1042 Fax: 800 762 6311

Email: questions@sekisui-dx.com



^{**}LOD for Influenza A/Michigan/45/15 (H1N1). Refer to IFU for additional LODs.







READ RESULTS AT 10 MINUTES