



QuickVue®


RSV



Test and Treat Today®!

- Results in 15 minutes, or less.
- All components necessary to perform test are contained in kit.
- Room temperature storage (24 months from date of manufacture).
- Built-in Positive and Negative controls confirm the test has performed properly.
- Positive and Negative external controls provided in the kit.
- Individually packaged for convenient single-use testing and disposal.

Nasopharyngeal aspirate and Nasal/Nasopharyngeal wash procedure:



Nasopharyngeal Aspirate:

Sensitivity: 99%

Specificity: 92%

Accuracy: 94%


Nasopharyngeal Wash:

Sensitivity: 83%

Specificity: 90%

Accuracy: 87%

Nasopharyngeal swab procedure:



Nasopharyngeal Swab:

Sensitivity: 92%

Specificity: 92%

Accuracy: 94%

The performance of any rapid RSV test is dependent on sample collection and handling and the adherence to the Package Insert. Refer to our website at quidel.com for additional performance claims.