

Cardinal Health™ Monoject™ Sharps Safety Product Selection Guide

2023-2024 vaccine administration

Products for your facility's unique vaccine delivery protocols, engineered to keep you safe


Follow this Product Selection Guide to quickly identify the products you'll need from medication prep through delivery and disposal, each designed to help keep you safe.



Determine the right products for your facility

Choose the path that aligns with your vaccine order type, CDC administration overviews and your facility's sharps protocol to discover the sharps safety products designed for your needs. Then locate corresponding product details and ordering information on page 2.

Vaccine order type



Draw

Order code	Description	Qty.
8881540111	Monoject™ SmarTip™ Vial Access Cannula	1,000/cs
8881540133	Monoject™ SmarTip™ 3 mL Vial Access Combination	800/cs
11811022	Monoject™ Blunt Fill Needle, 18 G x 1.5 in. (1.27 x 38 mm)	1,000/cs
1180100777	Monoject™ 1 mL Syringe Luer-Lock Tip	240/cs
1180300777	Monoject™ 3 mL Syringe Luer-Lock Tip	800/cs

Draw and inject

Order code	Description	Qty.
Magellan™ Hypodermic Safety Needle and Syringe Combination		
8881811310	1 mL, 23 G x 1 in.	50/box, 200/cs
8881811510	1 mL, 25 G x 1 in.	50/box, 200/cs
8881833810	3 mL, 18 G x 1 in.	50/box, 400/cs
8881833115	3 mL, 21 G x 1.5 in.	50/box, 400/cs
8881833210	3 mL, 22 G x 1 in.	50/box, 400/cs
8881833215	3 mL, 22 G x 1.5 in.	50/box, 400/cs
8881833310	3 mL, 23 G x 1 in.	50/box, 400/cs
8881833558	3 mL, 25 G x 5/8 in.	50/box, 400/cs
8881833510	3 mL, 25 G x 1 in.	50/box, 400/cs
Monoject™ Hypodermic Safety Needle and Syringe Combination		
11812558	1 mL, 25 G x 5/8 in.	100/box, 1,200/cs
11832110	3 mL, 21 G x 1 in.	100/box, 1,200/cs
11832115	3 mL, 21 G x 1.5 in.	100/box, 1,200/cs
11832215	3 mL, 22 G x 1.5 in.	100/box, 1,200/cs
11832310	3 mL, 23 G x 1 in.	100/box, 1,200/cs
11832315	3 mL, 23 G x 1.5 in.	100/box, 1,200/cs
11832558	3 mL, 25 G x 5/8 in.	100/box, 1,200/cs
11832510	3 mL, 25 G x 1 in.	100/box, 1,200/cs

Ordering tip:

Stock a variety of needle lengths

Because vaccines are commonly administered through intramuscular injections, the length of needle depends on a person's age/body size. Consider stocking a variety of needles in lengths of 5/8 in., 1 in. and 1.5 in. (22 G, 23 G or 25 G) to accommodate various ages and sizes.

Inject

Order code	Description	Qty.
Magellan™ Hypodermic Safety Needle (needle only)		
8881850110	21 G x 1 in.	50/box, 500/cs
8881850115	21 G x 1.5 in.	50/box, 500/cs
8881850210	22 G x 1 in.	50/box, 500/cs
8881850215	22 G x 1.5 in.	50/box, 500/cs
8881850358	23 G x 5/8 in.	50/box, 500/cs
8881850310	23 G x 1 in.	50/box, 500/cs
8881850558	25 G x 5/8 in.	50/box, 500/cs
8881850515	25 G x 1.5 in.	50/box, 500/cs
Monoject™ Hypodermic Safety Needle (needle only)		
1182110	21 G x 1 in.	100/box, 800/cs
1182115	21 G x 1.5 in.	100/box, 800/cs
1182215	22 G x 1.5 in.	100/box, 800/cs
1182310	23 G x 1 in.	100/box, 800/cs
1182315	23 G x 1.5 in.	100/box, 800/cs
1182558	25 G x 5/8 in.	100/box, 800/cs
1182510	25 G x 1 in.	100/box, 800/cs
1182515	25 G x 1.5 in.	100/box, 800/cs

Dispose

Order code	Description	Qty.
Monoject™ Sharps Disposal Containers and container hardware		
8506SA	In-room sharps container with SharpStar™ Lid, clear, 5 qt.	20/cs
8507SA	In-room sharps container with SharpStar™ Lid, transparent red, 5 qt.	20/cs
8515-	Tabletop holder for in-room sharps container, almond, 2 and 5 qt.	5/cs
85161H	Wall enclosure for in-room sharps container, 2 and 5 qt.	1/cs

To ensure your facility has the supplies it needs for vaccine administration, contact your sales representative, call **877.CARDINAL** or visit **cardinalhealth.com/sharpsafety/vaccineprepared**