INSTRUMENT REPROCESSING BROCHURE




Instrument Reprocessing Solutions Trusted by Thousands of Healthcare Facilities

Proper instrument reprocessing is vital to any facility's infection prevention program. Metrex offers precleaners, detergents, lubricants and associated products – an integrated system of products – that satisfy the exacting requirements for infection prevention to reduce the risk of cross-contamination.



The Instrument Reprocessing Workflow







Speeds Cleaning.

EmPower Foam dual enzymatic spray is excellent for dissolving bioburden immediately after surgery. Start in the operating room to prevent debris from drying on instruments.



Various Sizes.

Standard 1, 5, and 30-gallon sizes for manual and automated uses.



Saves Time.

EmPower dual-enzymatic cleaner saves time during the manual cleaning process by breaking down bioburden and keeping it off.



Compatibility.

An array of products work efficiently with a full range of washer/decontaminator/ sterilizer types.

All products featured will help you adhere to AAMI guidelines.



Metrex provides you with convenience and peace of mind.

Purchasing from multiple vendors costs extra time, money and space tied up in unnecessary inventory. One of our knowledgeable sales consultants can make it easy for you to get everything you need from one manufacturing source.

We'll help you...

- Get the best value for your dollar
- Ensure your washers are using the proper amounts of solutions
- Learn more about the importance of pre-cleaning through education tools

To find your local Metrex consultant, call 1.800.841.1428 or visit www.metrex.com.

EmPower™ Foam

Foaming Enzymatic Spray

A ready-to-use, foaming dual-enzymatic spray that starts the pre-cleaning process by breaking down blood, tissue, mucous and other protein-rich body fluids from instruments.

- Covers instruments, reducing splashes and spills
- Speeds up the cleaning process

Part No.Product Description10-4224EmPower Foam 24 ounce spray

Packaging 12 bottles/case

MetriSponge™

Pre-Moistened Enzymatic Sponge

MetriSponge is used to pre-clean critical or semi-critical medical devices prior to terminal sterilization or high level disinfection. This easy-to-use sponge, pre-moistened with MetriZyme dual-enzymatic detergent, helps ensure thorough pre-cleaning.

- Simplifies pre-cleaning of medical instruments
- Unique shape easily forms around cylindrical instruments
- Gentle on delicate instruments

Part No.	Product Description	Packaging
10-4025	MetriSponge - 25 sponges/box	4 boxes / case

EmPower™

Dual-Enzymatic Detergent

A dual-enzymatic detergent that works fast to remove blood, tissue, mucous and other protein-rich body fluids from instruments.

- · Contains two proteolytic enzymes for effective cleaning
- · Low-foaming formula is ideal for use in automated machines

Part No.	Product Description	Packaging
10-4102	EmPower – 2 ounces	48 bottles / case
10-4100	EmPower – 1 gallon	4 bottles / case
10-4150	EmPower – 5 gallon	1 drum
10-4130	EmPower – 30 gallon	1 drum
10-4400	EmPower Fragrance Free -1 gallon	4 bottles/case

MetriZyme™

Dual-Enzymatic Detergent

An extra-strength dual-enzymatic detergent that works fast to remove blood, tissue, mucous and other protein-rich body fluids from instruments.

· Low-foaming formula is ideal for use in automated machines

Part No.	Product Description	Packaging
10-4005	MetriZyme – 1 quart	4 bottles / case
10-4000	MetriZyme – 1 gallon	4 bottles / case

MetriClean[™]2

Multi-Purpose Instrument Cleaner/Polisher and Cart Washer Agent

An alkaline pH, low-foaming, concentrated cleaner for stainless steel instruments and cart washers.

- Removes blood and inorganic soils
- Removes stains and restores shine
- Increases instrument functionality

Part No.	Product Description	Packaging
10-8100	MetriClean2 – 1 gallon	4 bottles / case

MetriLube™

Instrument Lubricant Concentrate

An oil-based instrument lubricant that helps prevent spotting, staining and rusting of surgical instruments.

- Extends life of instruments
- Anticorrosive
- Nongreasy and silicone-free

Part No.	Product Description	Packaging
10-3425	MetriLube Spray – 24 ounce	12 bottles / case
10-3400	MetriLube – 1 gallon	4 bottles / case

MetriShine™

Descaler and Rust Remover

A two-in-one product that enhances the life and efficiency of instruments, automated washers and equipment.

- Removes water scale and mineral deposits
- Removes rust stains, spotting and corrosion from stainless steel instruments
- Saves money by restoring instruments

Part No.	Product Description	Packaging
10-9300	MetriShine – 1 gallon	4 bottles / case

Accessories

Part No.	Product Description	Packaging
46-1020	Dispensing Pump (quart)	each
46-1005	Dispensing Pump (1/2 gallon)	each
46-1000	Dispensing Pump (gallon)	each